

奇×助行器

Wheelchair x Rollator 2-in-1

1. 照顧者陪伴時...

可以推到平地讓長者自行運動 When going out with carer



輪椅 ←→ 助行器 /heelchair Rollator

迅速變身 鼓勵運動

Encourage exercise

2. 長者夫妻生活時...

個坐,一個推,一起度過夫妻生活 When travelling as a couple



輪椅 + 助行器 /heelchair Rollator Wheelchair

夫妻生活 輪流運動

Quality time together

3. 康復階段時...

可以一物兩用, 看護著不同康復階段 **During rehabilitation**



輪椅 → 助行器 Wheelchair

一物兩用 添加用途 Change usage as you recover





規格 Specifications	Rollz Motion ²	Rollz Motion ² Small
助行器總重量 Rollator wgt.	11.6 kg	11.6 kg
輪椅總重量 Wheelchair wgt.	15 kg	15 kg
負載 Max. Load	125 kg	125 kg
建議使用人士身高 User hgt.	1.60 - 1.90 m	1.50 - 1.85 m
座椅總高 Specifications	55 cm	50 cm
手推柄總高 Handle hgt.	87-98 cm	82-94 cm
助行器座椅尺寸 Rollator dim.	20 x 45 cm	20 x 45 cm
輪椅座椅尺寸 Wheelchair dim.	42 x 45 cm	42 x 45 cm
總閥 Frame width	67 cm	67 cm
總長 Frame length	63 cm	63 cm
摺合後總閣 Folded width	31 cm	31 cm
前輪尺寸 Front wheel dim.	8 inch	8 inch
後輪尺寸 Rear wheel dim.	12 inch	12 inch



澤 科 技 有 限 公 司 Delta Pyramax Co., Ltd.

香港九龍灣啟祥道22號開達大廈3樓A室 3/F Unit A, Kader Building, 22 Kai Cheung Road,

Kowloon Bay, Hong Kong. www.deltapyramax.com

Tel: (852) 2511 2118 Fax: (852) 2507 5078 Email: Sales@dpx.hk



takes you further

Product features Rollz

輪椅 × 助行器

— — 阿述罗 Wheelchair x Rollator 2-in

Award winning Dutch design rollator and wheelchair in one, for going out in style.

Stylish

安全可靠

Safe

後鼓式刹車,不怕路面濕滑

Innovative drum brakes, suitable to wet and dry conditions.

方便使用

Convenient

Rest on a comfortable cushioned seat or transform the rollator into a wheelchair.

易用輔助

Easy aid

Cross doorsteps and kerbs easily with the kerb assistant.



結構輕巧,前輪靈敏 Maneuverable

Light steering and easy motion with the maneuverable front wheels.

Walk tall and steadily between the ergonomic hand grips that are adjustable in height.

輕巧便携 Compact

容易摺合収納

The Rollz is easily foldable into a compact package that fits in almost any car.

Comfortable 出行及坐得更舒服

The large rear wheels with soft tyres and foam handles ensure a comfortable walk.













