

Asahi Industry produces a new urination care product



I always feel secure because I have Dandy Urinar.

★Supports independent urination★

# Dandy Urinar


Patent applied for / Design registration applied for

*A wearable urine collector that can be worn to sleep*

Reduces the burden of assistance required for night-time urination.



A proprietary idea that is easy and comfortable to use!



**1** The wavy shape of the inner receiver prevents sticking to the skin. Includes "barbs" which are cushioned and prevent backflow

**2** Medical-grade gentle texture made of soft silicone 2-layer design to maximize both fit and breathability



A large capacity urine tank (2.8 liters) Also works for nocturia, eliminating the need to handle frequent urination The cap is designed to be removable and attachable with one touch



Really easy! The surface fastener allows easy removal and attachment with no need for raising the hips Stays in place solidly with the specialized fixed garter  
※Comes equipped standard

**A great helper for those who provide care!**

It reduces the burden not only for the person with urinary troubles, but also the people who care for that person.

Now I can sleep soundly all night too!

I don't have to worry about inconveniencing anyone.

I don't have to call the nurse to urinate!



同梱セット



Accessories  
Set contents  
Receiver (urine reception implement) x 1  
Urine tube 1.5m  
Urine collection tank 2.8 l  
Cap  
1 specialized garter  
Backflow prevention membrane  
Absorbent sheets x 3  
Cleaning brush

※The specifications may be changed without notice for product improvements.

 **Asahi Industry Co., Ltd.**

1-11-22, Yonban, Atsuta-Ku, Nagoya, Aichi, Japan

TEL : +81-52-671-5193 FAX : +81-52-671-5196

E-mail : info\_e@asahi-sg.co.jp H.P. : http://asahi-packaging.com



**佳澤科技有限公司**

**Delta Pyramax Co., Ltd.**

香港九龍灣啟祥道22號開達大廈3樓A室  
3/F Unit A, Kader Building, 22 Kai Cheung Road, Kowloon Bay, Hong Kong.

Tel: (852) 2511 2118  
Fax: (852) 2507 5078  
Email: sales@dpx.hk